

Get full information about the ECG Dongle device
and the service "CardioCloud"

Buy ECG Dongle and supplies

Get service support

www.cardio-cloud.ru

8 499 130-98-92

Personal doctor:

Personal doctor ID:



Moscow, Varshavskoye highway, 125
+7 499 130-98-92

Moscow region, Dubna, Programmistov st., 4
+7 499 608-01-88

ECG Dongle

DIAGNOSTIC CARDIO COMPLEX

Cardio dongle

Mobile app

Cloud service



Professional guide to the world of your heart

1

Download the ECG Dongle mobile app



App Store


Google Play



2

Fix the electrodes on the body


Four reusable electrodes have adhesive surface and are easily placed on the body independently according to a colour indication.



3

Connect ECG Dongle to the mobile device

ECG Dongle lets you easily and accurately gets the information on the work of your heart at any time and in any place.



4

Take an examination

Go to the application and follow the prompts, record from 30 sec. to 10 min.



5

Send the data to "CardioCloud" or save the results to pdf format

In the "CardioCloud" service a cardiologist analyses the data and provides an opinion on the work of your heart and gives you recommendations.



CardioCloud

www.cardio-cloud.ru